

NEW!



CryOmega®

Disposable Cryosurgical Device

**Clinically effective
direct application
of liquid nitrous oxide
(-89°C) in a convenient
and economical device**




Premier® Medical Products
www.premusa.com

“I believe the CryOmega® device is appropriate for all Primary Care Physicians. It’s a great way to introduce cryosurgery for physicians who don’t realize how important full-time cryosurgery is for their practice.”

John E. Hocutt, Jr. M.D., FAAFP

It’s as easy as

1



2



3



Effective – Delivers a continuous, direct spray of nitrous oxide (-89°C) to treatment site. With a lower operating temperature compared to canister-based systems, CryOmega provides clinically effective freezing with pinpoint precision.

Convenient – Simple, hassle-free operation. Activate internal cartridge once and the device is ready for use. No need for additional accessories or replacement cartridges. 30 to 45 five-second sprays per device. Difficult lesions may require longer treatment.

Economical – Low-cost, effective cryosurgery in any clinical setting.

Disposable – Simply dispose of the entire unit when the internal cartridge is empty.

9034000 **CryOmega**

Disposable cryosurgical unit containing 16g nitrous oxide (N₂O) cartridge

CryOmega is approved for the following indications:

Verruca Vulgaris, Verruca Plantaris, Mollusca Contagiosa, Skin Tags, Lentigo, Actinic Keratosis, Verruca Plana, Seborrheic Keratosis and Genital Lesions

Please consult available medical literature to become fully informed on the benefits and risks of cryosurgical procedures.

	Before Treatment	Day 6	Day 10	Day 13
Treatment of Verruca Vulgaris <i>with one 13 second treatment</i>				

Available from:



Premier® Medical Products

1710 Romano Drive, Plymouth Meeting, PA 19462 U.S.A.
 610-239-6053 • 888 PREM USA (773-6872) • Fax: 610-239-6171
 medinfo@premusa.com • www.premusa.com