

Baxter

Welch Allyn

Spot

VISION SCREENER



PRODUCT SPECIFICATIONS

Dimensions	8 1/2 in x 6 3/4 in x 4 3/4 in (21.6 cm x 17.1 cm x 12.1 cm)
Weight	~2.55 lbs (1157 g)
Wireless Network	802.11 a/b/g/n
Operating Temperatures	50°F to 104°F (+10°C to +40°C)
Operating Humidity	30% to 75% relative humidity (non-condensing)
Storage/Shipping Temperatures	32°F to 122°F (0°C to +50°C)
Storage/Shipping Humidity	0% to 95% relative humidity (non-condensing)
Storage/Shipping Atmospheric Pressure	800 hPA to 1060 hPA
Expected Battery Life	2.5 years (under normal usage)

MEASUREMENTS

Spherical Equivalent	Range: -7.50D to +7.50D in 0.25D increments Accuracy: -3.50D to 3.50D, ± 0.50D; -7.50D to < -3.50D, ± 1.00D; > 3.50D to 7.50D, ± 1.00D
Cylindrical	Range: -3.00D to +3.00D in 0.25D increments Accuracy: -1.50D to 1.50D, ± 0.50D; -3.00D to < -1.50D, ± 1.00D; > 1.50D to 3.00D, ± 1.00D
Cylindrical Axis	Range: 1 to 180 degrees in 1 degree increments Accuracy: ± 10 degrees (for cylinder values > 0.5D)
Pupil Size	3.0 mm to 9.0 mm in 0.1 mm increments in ± cylinder* Accuracy: ± 0.4 mm
Pupil Distance	Range: 35 mm to 80 mm in 1 mm increments Accuracy: ± 1.5 mm
Gaze Measurement	0 to 20 degrees in increments of 1 degree nasal, temporal, inferior and superior Accuracy: ± 1.5 degrees
Time Per Measurement	Second(s)
Measuring Distance	3.3 ft ± 2 in (1 m ± 5 cm)
Fixation Target	Random visual pattern and audible sound
Measurement Principal	Optical reflex

For more information, please contact your local distributor or a Baxter sales representative at **1-800-535-6663**.

*3.0 mm available for 20-100 age range.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, Hillrom, Spot and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC158-230058 V1 08/23